

Аналитические колонки PolyRP

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: sxc@nt-rt.ru || сайт: <https://sepax.nt-rt.ru/>

Part	Product Name
------	--------------

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
2103 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	30		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	50		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	100		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	150		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
2125 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	250		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	30		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	50		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	100		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	150		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å,
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	250		

Part	Product Name
------	--------------

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
7830 7.8 x 300 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	7.8
Length (mm)	300		

266950- HPLC Column, Sepax, PolyRP-1000, 15 µm, 1000
21205 Å, 21.2 x 50 mm

Separation Mode	RP	Particle Size (µm)	15
Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	50		

266100- HPLC Column, Sepax, PolyRP-100, 15 µm, 100 Å,
21205 21.2 x 50 mm

Separation Mode	RP	Particle Size (µm)	15
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	50		

266100- HPLC Column, Sepax, PolyRP-100, 15 µm, 100 Å,
21230 21.2 x 300 mm

Separation Mode	RP	Particle Size (µm)	15
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	300		

266950- HPLC Column, Sepax, PolyRP-1000, 15 µm, 1000
21230 Å, 21.2 x 300 mm

Separation Mode	RP	Particle Size (µm)	15
Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	300		

261950- Conventional Column, Sepax, PolyRP-1000, 10
50025 µm, 1000 Å 50.0 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	1000	ID (mm)	50
Length (mm)	250		

261950F- FLOW Axial Column, Sepax, PolyRP-1000, 10 µm,
50025 1000 Å 50.0 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	50
Length (mm)	250		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
21230 21.2 x 300 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	300		

261135- HPLC Column, Sepax, PolyRP-Mix, 10 µm, Mixed
7830 Pore, 7.8 x 300 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	Mix	ID (mm)	7.8
Length (mm)	300		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
21230 21.2 x 300 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	21.2
Length (mm)	300		

Part	Product Name
------	--------------

261200- HPLC Column, Sepax, PolyRP-100/300, 10 µm,
 4615 Mixed Bed 100 Å and 300 Å, 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100 and 30	ID (mm)	4.6
Length (mm)	150		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å,
 4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	100		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
 7830 7.8 x 300 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	7.8
Length (mm)	300		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
 7805 7.8 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	7.8
Length (mm)	50		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
 7805 7.8 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	7.8
Length (mm)	50		

260300- HPLC Guard Column Plus Holder, Sepax, PolyRP-
 21201C 300, 5 µm, 300 Å 21.2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	10		

262005- HPLC Guard Column Plus Holder, Sepax, PolyRP-
21201C NP5, 5 µm, NP 21.2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	10		

262010- HPLC Guard Column Plus Holder, Sepax, PolyRP-
21201C NP10, 10 µm, NP 21.2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	10		

261300- HPLC Guard Column Plus Holder, Sepax, PolyRP-
21201C 300, 10 µm, 300 Å 21.2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	10		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
50025 50 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	50
Length (mm)	250		

Part	Product Name
------	--------------

261500- PolyRP-500, 10 µm 500 Å 50 x 300 mm LL
 50030 (packed with 150 mm polyRP)

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	50
Length (mm)	300		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
 7825 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	7.8
Length (mm)	250		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 2.1
 2105 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	50		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
 4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	250		

262005- HPLC Column, Sepax, PolyRP-NP5, Guard
 2001C cartridge with holder, 5 µm, NP 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
 4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	150		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 2 x
2001 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	50		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 10
10025 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	10
Length (mm)	250		

262010- HPLC Column, Sepax, PolyRP-NP10, Guard
4001C cartridge with holder, 10 µm, NP 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		

Part	Product Name
------	--------------

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP 4
4001 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	250		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	150		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	100		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	150		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	50		

262010- HPLC Column, Sepax, PolyRP-NP10, Guard
2001C cartridge with holder, 10 µm, NP 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP 2
2001 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	10
Length (mm)	250		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 7.8
7825 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	7.8
Length (mm)	250		

Part	Product Name
------	--------------

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 4.6
4625 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	250		

262003- HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 4.6
4625 x 250 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	250		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 4.6
4615 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	150		

262003- HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 4.6
4615 x 150 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	150		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 4.6
4605 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	50		

262003- HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 4.6
4610 x 100 mm

Separation Mode	RP	Particle Size (µm)	3
------------------------	----	---------------------------	---

Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	100		
262005-4001C	HPLC Column, Sepax, PolyRP-NP5, Guard cartridge with holder, 5 µm, NP 4 x 10 mm		
Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		
262003-4605	HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 4.6 x 50 mm		
Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	50		
262005-4001	HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 4 x 10 mm		
Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		
262003-4603	HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 4.6 x 30 mm		
Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	4.5
Length (mm)	30		

Part	Product Name
------	--------------

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	250		

262003- HPLC Column, Sepax, PolyRP-NP3, Guard
4001C cartridge with holder, 3 µm, NP 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	150		

262003- HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 4 x
4001 10 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	NP	ID (mm)	21.2
Length (mm)	100		

262003- HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 2.1
2115 x 150 mm

Separation Mode	RP	Particle Size (µm)	3
------------------------	----	---------------------------	---

Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	150		

262005- HPLC Column, Sepax, PolyRP-NP5, 5 µm, NP 2.1
2115 x 150 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	150		

262003- HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 2.1
2105 x 50 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	50		

262003- HPLC Column, Sepax, PolyRP-NP3, Guard
2001C cartridge with holder, 3 µm, NP 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

262003- HPLC Column, Sepax, PolyRP-NP3, 3 µm, NP 2 x
2001 10 mm

Separation Mode	RP	Particle Size (µm)	3
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

Part	Product Name
------	--------------

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
7825 Å 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	7.8
Length (mm)	250		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
7815 Å 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	7.8
Length (mm)	150		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	250		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
4625 Å 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	250		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	150		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
4615 Å 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	150		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	100		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
4610 Å 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	100		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	50		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
4605 Å 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	50		

Part	Product Name
------	--------------

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	30		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
4603 Å 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	30		

262002- HPLC Column, Sepax, PolyRP-NP1.7, Guard
4001C cartridge with holder, 1.7 µm, NP 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		

261950- HPLC Column, Sepax, PolyRP-1000, Guard
4001C cartridge with holder , 10 µm, 1000 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	4
Length (mm)	10		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
4001 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	4
Length (mm)	10		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
4001 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	1000	ID (mm)	4
Length (mm)	10		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
30025 Å 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	30
Length (mm)	250		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	150		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
30015 Å 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	30
Length (mm)	150		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	50		

Part	Product Name
------	--------------

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
2125 Å 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	250		

262002- HPLC Column, Sepax, PolyRP-NP1.7, Guard
2001C cartridge with holder, 1.7 µm, NP 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
21225 Å 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	250		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
2001 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	2
Length (mm)	10		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
21215 Å 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	150		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
21210 Å 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	100		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
2115 Å 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	150		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
2110 Å 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	100		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
2105 Å 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	50		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
2103 Å 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	30		

Part	Product Name
------	--------------

261950- HPLC Column, Sepax, PolyRP-1000, Guard
2001C cartridge with holder , 10 µm, 1000 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	2
Length (mm)	10		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
2001 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	2
Length (mm)	10		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
10025 Å 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	10
Length (mm)	250		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
10015 Å 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	10
Length (mm)	150		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
10010 Å 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	1000	ID (mm)	10
Length (mm)	100		

261950- HPLC Column, Sepax, PolyRP-1000, 10 µm, 1000
10005 Å 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	1000	ID (mm)	10
Length (mm)	50		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
7825 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	7.8
Length (mm)	250		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
7815 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	7.8
Length (mm)	150		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	250		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	150		

Part	Product Name
------	--------------

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	100		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	50		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	30		

261500- HPLC Column, Sepax, PolyRP-500, Guard
4001C cartridge with holder , 10 µm, 500 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	4
Length (mm)	10		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
4001 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	4
Length (mm)	10		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
4001 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	300	ID (mm)	4
Length (mm)	10		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
30025 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	30
Length (mm)	250		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
30025 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	30
Length (mm)	250		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
30015 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	30
Length (mm)	150		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
30015 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	30
Length (mm)	150		

Part	Product Name
------	--------------

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
2125 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	250		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
2125 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	250		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	21.2
Length (mm)	250		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	250		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	21.2
Length (mm)	150		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	150		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	21.2
Length (mm)	100		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	100		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	150		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	150		

Part	Product Name
------	--------------

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	100		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	100		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	50		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	50		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
2103 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	30		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
2103 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	30		

261500- HPLC Column, Sepax, PolyRP-500, Guard
2001C cartridge with holder , 10 µm, 500 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	2
Length (mm)	10		

261300- HPLC Column, Sepax, PolyRP-300, Guard
2001C cartridge with holder , 10 µm, 300 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	2
Length (mm)	10		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
2001 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	2
Length (mm)	10		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
2001 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	2
Length (mm)	10		

Part	Product Name
------	--------------

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	10
Length (mm)	250		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	10
Length (mm)	250		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	10
Length (mm)	150		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	10
Length (mm)	150		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
10010 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	10
Length (mm)	100		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
10010 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	300	ID (mm)	10
Length (mm)	100		

261500- HPLC Column, Sepax, PolyRP-500, 10 µm, 500 Å
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	500	ID (mm)	10
Length (mm)	50		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	10
Length (mm)	50		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
7825 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	7.8
Length (mm)	250		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
7825 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	7.8
Length (mm)	250		

Part	Product Name
------	--------------

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
7815 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	7.8
Length (mm)	150		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
7815 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	7.8
Length (mm)	150		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	250		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	250		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	150		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	150		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	100		

261200- HPLC Column, Sepax, PolyRP-100/300, 10 µm,
4610 Mixed Bed 100 Å and 300 Å, 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100 and 30	ID (mm)	4.6
Length (mm)	100		

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	50		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	50		

Part	Product Name
------	--------------

261300- HPLC Column, Sepax, PolyRP-300, 10 µm, 300 Å
 4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	30		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
 4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	30		

261300- HPLC Column, Sepax, PolyRP-300, Guard
 4001C cartridge with holder , 10 µm, 300 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	300	ID (mm)	4
Length (mm)	10		

261100- HPLC Column, Sepax, PolyRP-100, Guard
 4001C cartridge with holder, 10 µm, 100 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	4
Length (mm)	10		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
 4001 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	4
Length (mm)	10		

260110100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
 30025 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
------------------------	----	---------------------------	----

Pore Size (Å)	100	ID (mm)	30
Length (mm)	250		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
30015 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	30
Length (mm)	150		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
2125 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	250		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	250		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	150		

Part	Product Name
------	--------------

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	100		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	150		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	100		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	50		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
2103 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	2.1
Length (mm)	30		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
7825 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	1000	ID (mm)	7.8
Length (mm)	250		

261100- HPLC Column, Sepax, PolyRP-100, Guard
2001C cartridge with holder, 10 µm, 100 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	2
Length (mm)	10		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
7815 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	7.8
Length (mm)	150		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
2001 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	2
Length (mm)	10		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	250		

Part	Product Name
------	--------------

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	10
Length (mm)	250		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	150		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	10
Length (mm)	150		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	100		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
10010 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	10
Length (mm)	100		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	50		

261100- HPLC Column, Sepax, PolyRP-100, 10 µm, 100 Å
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	10
Length (mm)	50		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	4.6
Length (mm)	30		

260950- HPLC Column, Sepax, PolyRP-1000, Guard
4001C cartridge with holder , 5 µm, 1000 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	4
Length (mm)	10		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
4001 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	4
Length (mm)	10		

Part Product Name

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
30025 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	30
Length (mm)	250		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
30015 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	30
Length (mm)	150		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
2125 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	250		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	250		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	150		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	1000	ID (mm)	21.2
Length (mm)	100		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	150		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	100		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	50		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
2103 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	2.1
Length (mm)	30		

Part	Product Name
------	--------------

260950- HPLC Column, Sepax, PolyRP-1000, Guard
2001C cartridge with holder , 5 µm, 1000 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	2
Length (mm)	10		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
2001 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	2
Length (mm)	10		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	10
Length (mm)	250		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	10
Length (mm)	150		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
10010 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	10
Length (mm)	100		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	1000	ID (mm)	10
Length (mm)	50		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
7815 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	7.8
Length (mm)	150		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
7815 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	7.8
Length (mm)	150		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	250		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
4625 4.6 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	250		

Part	Product Name
------	--------------

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	150		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
4615 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	150		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	100		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	100		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	50		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	50		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	4.6
Length (mm)	30		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
4603 4.6 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	4.6
Length (mm)	30		

260500- HPLC Column, Sepax, PolyRP-500, Guard
4001C cartridge with holder , 5 µm, 500 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	4
Length (mm)	10		

260300- HPLC Column, Sepax, PolyRP-300, Guard
4001C cartridge with holder , 5 µm, 300 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	4
Length (mm)	10		

Part	Product Name
------	--------------

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å 4
4001 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	4
Length (mm)	10		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å 4
4001 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	4
Length (mm)	10		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
30025 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	30
Length (mm)	250		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
30025 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	30
Length (mm)	250		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
30015 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	30
Length (mm)	150		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
30015 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	300	ID (mm)	30
Length (mm)	150		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
2125 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	250		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
2125 2.1 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	250		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	21.2
Length (mm)	250		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	250		

Part	Product Name
------	--------------

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	21.2
Length (mm)	150		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	150		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	21.2
Length (mm)	100		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	21.2
Length (mm)	100		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	150		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
2115 2.1 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	150		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	100		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	100		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	50		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
2105 2.1 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	50		

Part	Product Name
------	--------------

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
 2103 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	2.1
Length (mm)	300		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
 2103 2.1 x 30 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	2.1
Length (mm)	30		

260500- HPLC Column, Sepax, PolyRP-500, Guard
 2001C cartridge with holder , 5 µm, 500 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	2
Length (mm)	10		

260300- HPLC Column, Sepax, PolyRP-300, Guard
 2001C cartridge with holder , 5 µm, 300 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	2
Length (mm)	10		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å 2
 2001 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	2
Length (mm)	10		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å 2
 2001 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	300	ID (mm)	2
Length (mm)	10		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	10
Length (mm)	250		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	10
Length (mm)	250		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	10
Length (mm)	150		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	10
Length (mm)	100		

Part	Product Name
------	--------------

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
10010 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	10
Length (mm)	100		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
10010 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	10
Length (mm)	100		

260500- HPLC Column, Sepax, PolyRP-500, 5 µm, 500 Å
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	10
Length (mm)	50		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	10
Length (mm)	50		

260300- HPLC Column, Sepax, PolyRP-300, 5 µm, 300 Å
7825 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	300	ID (mm)	7.8
Length (mm)	250		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
7825 7.8 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	100	ID (mm)	7.8
Length (mm)	250		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
7815 7.8 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	7.8
Length (mm)	150		

260100- HPLC Column, Sepax, PolyRP-100, Guard
4001C cartridge with holder, 5 µm, 100 Å 4 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	4
Length (mm)	10		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å 4
4001 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	4
Length (mm)	10		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
30025 30 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	30
Length (mm)	250		

Part	Product Name
------	--------------

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
30015 30 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	30
Length (mm)	150		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
21225 21.2 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	250		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
21215 21.2 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	150		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
21210 21.2 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	100		

260100- HPLC Column, Sepax, PolyRP-100, Guard
2001C cartridge with holder, 5 µm, 100 Å 2 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	2
Length (mm)	10		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å 2
2001 x 10 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

Pore Size (Å)	100	ID (mm)	2
Length (mm)	10		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	10
Length (mm)	250		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	10
Length (mm)	150		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
10010 10 x 100 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	10
Length (mm)	100		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	10
Length (mm)	50		

Part	Product Name
------	--------------

263100- HPLC Column, Sepax, PolyRP-100, 30 µm, 100 Å,
4605 4.6 x 50 mm

Separation Mode	RP	Particle Size (µm)	30
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	50		

263100- HPLC Column, Sepax, PolyRP-100, 30 µm, 100 Å,
21205 21.2 x 50 mm

Separation Mode	RP	Particle Size (µm)	30
Pore Size (Å)	100	ID (mm)	21.2
Length (mm)	50		

263100- HPLC Column, Sepax, PolyRP-100, 30 µm, 100 Å,
10005 10 x 50 mm

Separation Mode	RP	Particle Size (µm)	30
Pore Size (Å)	100	ID (mm)	10
Length (mm)	50		

260100- HPLC Column, Sepax, PolyRP-100, 5 µm, 100 Å
7830 7.8 x 300 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	7.8
Length (mm)	300		

260950- HPLC Column, Sepax, PolyRP-1000, 5 µm, 1000 Å
7830 7.8 x 300 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	1000	ID (mm)	7.8
Length (mm)	300		

260115300- HPLC Column, Sepax, PolyRP-300, 15 µm, 300 Å
10015 10 x 150 mm

Separation Mode	RP	Particle Size (µm)	15
------------------------	----	---------------------------	----

Pore Size (Å)	300	ID (mm)	10
Length (mm)	150		

260115300- HPLC Column, Sepax, PolyRP-300, 15 µm, 300 Å
10025 10 x 250 mm

Separation Mode	RP	Particle Size (µm)	15
Pore Size (Å)	300	ID (mm)	10
Length (mm)	250		

262010- HPLC Column, Sepax, PolyRP-NP10, 10 µm, NP
4610 4.6 x 100 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	NP	ID (mm)	4.6
Length (mm)	10		

261100- FLOW Axial Column, Sepax, PolyRP-100, 10 µm,
50025 100 Å 50.0 x 250 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	50
Length (mm)	250		

266300- HPLC Column, Sepax, PolyRP-300, 15 µm, 300 Å,
30025 30.0 x 250 mm

Separation Mode	RP	Particle Size (µm)	15
Pore Size (Å)	300	ID (mm)	30
Length (mm)	250		

Part	Product Name
------	--------------

266500- HPLC Column, Sepax, PolyRP-500, 15 µm, 500 Å,
30025 30.0 x 250 mm

Separation Mode	RP	Particle Size (µm)	15
Pore Size (Å)	500	ID (mm)	30
Length (mm)	250		

260500- HPLC Column, Sepax, PolyRP, 5 µm, 500 Å 7.8 x
7830 300 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	500	ID (mm)	7.8
Length (mm)	300		

260100P- HPLC Column, Sepax, PolyRP-100, PEEK tubing,
4615 5 µm, 100 Å 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	5
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	150		

261100P- HPLC Column, Sepax, PolyRP-100, PEEK tubing,
4615 10 µm, 100 Å 4.6 x 150 mm

Separation Mode	RP	Particle Size (µm)	10
Pore Size (Å)	100	ID (mm)	4.6
Length (mm)	150		

262002- HPLC Column, Sepax, PolyRP-NP1.7, 1.7 µm, NP
2110 2.1 x 100 mm

Separation Mode	RP	Particle Size (µm)	1.7
Pore Size (Å)	NP	ID (mm)	2.1
Length (mm)	100		

260300- HPLC Column, Sepax, PolyRP, 5 µm, 300 Å 7.8 x
7830 300 mm

Separation Mode	RP	Particle Size (µm)	5
------------------------	----	---------------------------	---

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: sxc@nt-rt.ru || сайт: <https://sepax.nt-rt.ru/>